
 Poids maximum utilisateur/Max. user weight/Maximumgewicht gebruiker/Peso máximo utilizador/Peso max. Utente/Peso máximo do utilizador/Maximales Benutzergewicht  Poids du produit/Product's weight/Productgewicht/Peso del producto/Peso del prodotto/Peso do produto/Gewicht des Produktes

Référence P / Reference P / Referentie P / Código P / Codice P / Referência P / Referenz P <span> </span> :	047685 047686
Fonction Function Functie Función Funzione Função Funktion	Permet de prendre une douche en position assise. Allows to take a shower in a seated position. Zo kan u zittend een douche nemen. Permite tomar una ducha en posición sentada. Permette di prendere una doccia in posizione seduta. Permite tomar um duche sentado. Ermöglicht es, im Stzen zu duschen.
Descriptif Description Beschrijving Descripción Descrizione Função Beschreibung	 <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>150 kg</p> </div> <div style="text-align: center;">  <p>9 kg</p> </div> </div>
Numéro notice / Instruction sheet no / Fichenummer Número instrucciones / Numero manuale / Número do manual / Gebruiksaanwijzingsnummer	S00001833
Date de création / Creation's date / Aanmaakdatum / Fecha de creación Data di creazione / Data de criação / Erstellungsdatum	13/01/2011
Révision n° / Revision no. / Revisie nr. / Revisión n° / Revisione n° / Revisão n° / Revisionsnr.	02
Date / Date / Datum / Fecha / Data / Data / Datum	23/04/2013

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**GB  Do not install this equipment before reading the following instructions:**  **Fitting instructions.** We advise you to have your chair installed by a professional. He should first make sure that the entire surface of the wall to be covered by the chair is flat (no protruding plaster or tiles), and that the structure of the wall is suitable and, above all, does not display any weak spots. This will make it possible to define the best place to install the chair. We do not recommend Ergosoft chairs to be mounted on walls made of plasterboard, brick tiles or plaster tiles. The chair can only be installed on this type of wall if it can be fixed to a solid part or if specific reinforcements are used. This product is delivered ready-mounted. To install the chair properly, make sure the parts outlined in the assembly diagrams are removed. Before drilling into the wall, physically hold the product (in the lowered position) against the wall so that the stand is PERPENDICULAR to the floor, then mark the position of the holes to be drilled. Plugs have to be selected according to the type of wall. The 6 fixings points have to be necessarily respected.  **Care and maintenance.** IMPORTANT: In order to guarantee the lasting quality of this product and its safety for the user, we recommend that when cleaning it you should use mild soap solutions. **CAUTION: DO NOT USE** cleaning P = PELLET ASC

products made of: etones, hydrocarbons, acids, bases, esters and ethers. The disinfection has to be made according to the recommendations of use of the various products of disinfection. After use of a cleaning product, make a rinsing

in the clear water.  **Safety warnings.** Ergosoft chairs are designed to vertically bear maximum static loads of 150 kg. Each armrest can bear a load of 100 Kg. If these maximum loads are exceeded, the chair and armrests may be damaged. Ergosoft chairs should only be used as described in the manual. The materials used to manufacture these chairs allow flexibility when the chairs are used correctly. It is normal that the chair gives slightly in relation to its vertical position when it is bearing a load. The shower seat adds security to the user if it is fitted according to instructions. After fitting and before using, make sure that the fixing screws are tightened securely. The shower seat is designed to be used in seated position. To guarantee the user’s safety, it is not advised to rise up on the seat. When the seat is in a raised position, DO NOT USE it as a grab bar. To avoid any pinching risks, do not put fingers in the mobile parts. The product does not contain any hazardous substances. The manufacturer shall have no liability for any injury caused by improper installation or use of this shower seat.  At the end of its life cycle, the user must have the product recycled by a certified facility in the country of use.

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Fabricant : GODONNIER SAS